



Shanghai Joy Crown Industry Co., Ltd.

Add: Rm.1105 No.55-56 Lane 21 Pu San Road, Pu Dong New Area,

Shanghai 200125, P.R. China

Tel: +86-21-58751019 Fax: +86-21-50907607 Mobile/Whatsapp: 008613818337572

Email: sales@joycrown.com <https://www.medicalsupplyment.com>

Cheesecloth Gauze Specifications

Catalogue No. /Reference No. /Products Code	JC01015
Material	100% nature cotton gauze, unbleached and white, food grade, high absorbency, 100% safe and healthy
Cotton Yarn	40's cotton yarn
Gauze Mesh	19x8, 19x15, 20x12, 20x24, 28x16, etc
Size	36"x1yard/4ply/bag, 36"x2yds/4ply/bag, 36"x3yds/4ply/bag, 36"x4yds/4ply/bag, 36"x5yds/4ply/bag, 36"x6yds/4ply/bag, 36"x500yds/4ply/bag, 36"x1000yds/4ply/bag
Ply	4 Ply
Color	Unbleached and White
Sterile or Not	No Sterile
Rules & Classification	Rule 1, Class I





Shanghai Joy Crown Industry Co., Ltd.

Add: Rm.1105 No.55-56 Lane 21 Pu San Road, Pu Dong New Area,

Shanghai 200125, P.R. China

Tel: +86-21-58751019 Fax: +86-21-50907607 Mobile/Whatsapp: 008613818337572

Email: sales@joycrown.com <https://www.medicalsupplyment.com>

Cheesecloth Gauze Specifications

Reference No/ Product Code	Description	Size	Unit Package	Outside/Carton Package
JC010151	Cheese Cloth, 40's, 19x8 or other gauze meshes, different sizes, 1pc/PE bag with cardboard, with OEM printing on bag or label, 288/144 pcs/carton or 2 rolls/carton.	36" x 1yds-4ply	1 PC	288 PCS
JC010152		36" x 2yds-4ply	1 PC	288 PCS
JC010153		36" x 3yds-4ply	1 PC	144 PCS
JC010154		36" x 4yds-4ply	1 PC	144 PCS
JC010155		36" x 5yds-4ply	1 PC	144 PCS
JC010156		36" x 6yds-4ply	1 PC	144 PCS
JC010157		36" x 500yds-4ply	1 Roll	2 Rolls
JC010158		36" x 1,000yds-4ply	1 Roll	2 Rolls